Anti-Rabies Virus Mouse Monoclonal Antibody Clone 1112-1C8-3A5

Catalog number RDB.1112-1p500

Description

Anti-Rabies Virus Mouse Monoclonal Antibody
Species: Mouse monoclonal
Isotype: IgG1κ
Specificity: glycoprotein of ERA (Evelyn Rokitniki Abelseth) strain Rabies virus
Clone: 1112-1C8-3A5
Format: Unconjugated
Material size: 500 µg

Specifications

Source: *In vitro* production,
Purification: Protein A
Formulation: lyophilized
Purity : > 95 %
Qty: 500 µg

Storage instructions:
Pulse the vial before resuspension the antibody in ultrapure sterile water.
Store at +4°C for short term (1-2 weeks). Store at -20°C for long term.
Avoid repeated freeze / thaw cycles.

Tested applications: ELISA (qualified for sandwich ELISA with anti Rabies monoclonal antibody clone D1-25)
Optimal dilutions/concentrations should be determined by the end user.

Related products

Code	Product
RDB.1112-1p100	Anti-Rabies Virus, Mouse Monoclonal Antibody, Clone 1112-1C8-3A5, purified size 100 µg
RDB.D1-25p100	Anti-Rabies Virus, Mouse Monoclonal Antibody, Clone D1-25, purified 100 µg
RDB.D1-25b100	Anti-Rabies Virus, Mouse Monoclonal Antibody, Clone D1-25, biotinylated size 100 µg
RDB.D1-25p500	Anti-Rabies Virus, Mouse Monoclonal Antibody, Clone D1-25, purified 500 µg
RDB.D1-25b500	Anti-Rabies Virus, Mouse Monoclonal Antibody, Clone D1-25, biotinylated size 500 µg

Bibliography

Chabaud-Riou M, Moreno N, Guinchard F, Nicolai M-C, Niogret-Siohan E, Sève N, Manin C, Guinet-Morlot F & Riou P, Biologicals (2017), 46, 124-129
Dietzschold B, Rupprecht CE, Tollis M, Lafon M, Mattei J, Wiktor TJ, Koprowski H. Rev Infect Dis. 1988 Nov-Dec;10 Suppl 4:S785-98.
Lafon M, Perrin P, Versmisse & Sureau P, Journal of Biological Standardization (1985), 13, 295-301
Chabaud-Riou M, Moreno N, Guinchard F, Nicolai M-C, Niogret-Siohan E, Sève N, Manin C, Guinet-Morlot F & Riou P, Biologicals (2017), 46, 124-129
Müller T, Dietzschold B, Ertl H, Fooks AR, Freuling C, Fehlner-Gardiner C, et al. PLoS Negl Trop Dis 2009 Nov 3;3(11):e542.
Schumacher CL, Dietzschold B, Ertl HC, Niu HS, Rupprecht CE, Koprowski H. J Clin Invest 1989 Sep;84(3):971e5.